



Dino-Lite Medical



Hand Held Digital Microscope



The Dino-Lite is a compact digital microscope with USB - PC connectivity. Magnification ranges from 10X to 200X (adjustable single lens) to 500 X. LED lights are built-in with white, **polarized**, ultraviolet, infrared or switchable combinations depending on the model. Convenient and easy-to-use software is compatible with *Windows and Apple*. (Continuous magnification without changing lenses)
 Digital hand held Dermatoscopes have a closed and polarizing cap (through a disposable contact cap)
 Dino-Lite is ideal for recording detailed examinations of the skin with accurate measurement & calibration functions.

Majac Medical Products has been appointed the sole representative in Australia and New Zealand for the DinoLite Medical range of scopes. All products carry the CE Mark and are TGA listed.

DERMATOSCOPES

DERMATOSCOPE

Alloy Housing - **POLARIZED** - Resolution: **5M pixels**
 ON/OFF controlled by software – Micro touch snapshot
 Accurate measurement & calibration capability
 Capture a photograph, video or time-lapsed video with a button's click . Interface: USB (PC) 2.0
 Frame rate: up to 30fps, 15fps at 1,3M. Illumination: 8 white LED
 "AM35B" stand included Accessories included: P2 - Contact cap



Code: FEDMEDM700PAluminium



DERMATOSCOPE

Aluminium Alloy Housing Resolution: **5M pixels**
 Illumination: built-in 8 white LED
 ON/OFF controlled by software – Micro touch snapshot
 Accurate measurement & calibration capability
 Capture a photograph, video or time-lapsed video with a button's click. Interface: USB (PC) 2.0
 Frame rate: up to 30fps, 15fps at 1,3 M
 "AM35B" stand included. Accessories included: P1 - Contact cap



Code: FEDMEDM700



DERMATOSCOPE

Resolution: **1,3M pixels**
 UV Led on/off + WHITE LED on/off
 Medical Device Class1with Measurement Function
POLARIZED on/off Illumination: built - in 8 white LED
 ON/OFF controlled by software - Micro touch snapshot
 Accurate measurement & calibration capability
 Capture a photograph, video or time - lapsed video with a button's click. Interface: USB (PC) 2.0
 Frame rate: up to 30fps, 15fps at 1,3M - Interface: USB (PC) 2.0

POLARIZED on/off

Gode: AM413ZT



DERMATOSCOPE - As above with Aluminium Alloy Housing **POLARIZED on/off** Code:AM413MZT



DERMATOSCOPE

Resolution: 1,3M pixels Illumination: built in 8 white LED
ON/OFF controlled by software - Micro touch snapshot

Accurate measurement & calibration capability - Micro touch snapshot

Capture a photograph, video or time - lapsed video with a button's click

Frame rate: up to 30fps, 15fps at 1,3M - Interface: USB (PC) 2.0

DinoCapture Software included

"AMITA" stand included. Accessories included: P1 - Contact cap

Code: AM413T



DERMATOSCOPE - As above with Aluminium Alloy Housing code: AM413MT

DERMATOSCOPE FOR AESTHETIC USE

Resolution: 1,3M pixels

Illumination: built - in 4 WHITE LED

Illumination: built - in 4 UV LED

ON/OFF controlled by software - Micro touch snapshot

Excitation - 400Nm - Emission Filter Cut-off Wavelength - 425Nm

Accurate measurement & calibration capability

Capture a photograph, video or time - lapsed video

with a button's click. Interface: USB (PC) 2.0

Frame rate: up to 30fps, 15fps at 1,3M - Interface: USB (PC) 2.0

DinoCapture Software included

"AMITA" stand included. Accessories included: P1 - Contact cap

Code: AM413TFVW



DERMATOSCOPE AESTHETIC As above with Aluminium Alloy Housing Code: AM413MTFVW

DERMATOSCOPE MOBILE SOLUTIONS

The AMK412Z-C200 (Polarized) or AMK-412-C200 (not Polarized) are for mobile applications and has been designed as a dedicated combination of a Dino-Lite and a Portable Digital Recorder. The recorder has an 8.9cm LCD panel. Pictures and video can be recorded onto a SD card. Dedicated Dino-Lite models with Polarizer, Long Working Distance or Ultraviolet lighting are available.

Code: AMK - 412Z - C200

Resolution 640*480 pixels (VGA)

POLARIZED - OPTIONAL

Adjustable magnification: 10x ~ 50x, 200x

(Continuous magnification without switching lenses)

POLARIZED Built in 8 white - light LED illumination

Interface: AV, NTSC or PAL (TV)

AC plug in or 3.7 volt Rechargeable battery AC110/240V

Monitor LCD 3.5" - Built - in microphone and speaker

SD Flash Memory (not included)

Microscope Dimensions: 10cm (H) *3.2cm (D) Weight: 90g

LCD Dimensions :13cm*8cm*2,5cm - LCD Weight: 190g

Accurate Measurement Function

Capture a photograph, video or time - lapsed video with a button's click Accessories included: P2 - Contact cap





ACCESSORIES

Contact Cap (Small) for use in dermatology with the following Dermatoscopes

Dino-Lite Medical FEDMEDM700

Dino-Lite Medical AM413T - Dino-Lite Medical AM413TL

Dino-Lite Medical AM413TLFVW - Dino-Lite Medical AM413FIT

Dino-Lite Medical AM411T - Dino-Lite Medical AM413FVT

Dino-Lite Medical AM413FVT2 - Dino-Lite Medical AM413TFVW

Code: P1



Contact Cap (Large) for use in dermatology with the following Dermatoscopes

Dino-Lite Medical FEDMEDM700P **with Polarization**

Dino-Lite Medical AM413ZT **with Polarization**

Dino-Lite Medical AMK-412Z-C200 **with Polarization**

Code: P2



Material: Polycarbonate for medical use - Good hydrolysis resistance

Disinfection and Sterilization: Autoclave or ethelene oxide gas

TONGUE DEPRESSOR ATTACHMENT Code: MSAA501

The MSAA501 attachment allows you to attach a disposable tongue depressor onto a Dino - Lite. With this accessory you can inspect and document the larynx and do other oral inspections. Simply attach the MSAA501 to the Dino - Lite and then insert a tongue depressor and you are ready to examine. Compatible with all "AM" series *without Polarization*

FOOT CONTROL Code: SW-F1

The SW-F1 allows you to take a picture with a simple step from your foot using our Dino Lite or Dino Eye digital microscopy solutions. The foot pedal is compatible with both DinoXcope for Macintosh computers and DinoCapture 2.0 for PC, making it a single accessory for all users (Not compatible with DinoScope and prior version of DinoCapture).

There is no software CD required, simply plug in the device into your USB port and you can use it in conjunction with our Dino-Lite or Dino-Eye digital.

DESK TOP STANDS Code: AM 35B

The AM35B is an economical and rigid vertical Desk top stand that was designed to be used with the Dino Lite handheld, digital microscope cameras. It provides up to 360 degrees of horizontal axis rotational perspective and as much as 20cm of vertical working distance from the object.



Code: AM 35B



"AMITA" stand with flexible arm



Hand Held Digital Microscope



DERMATOSCOPY


AM413ZT

AM413ZT - X70



Skin check with Dinocapture

Dinocapture



AM413ZT with seal cap – 200X

AM413ZT with Polarize

