

LIBERTY 7E-D MULTI-VOLTAGE PORTABLE SUCTION PUMP

Liberty[®]



Features:

- Multi-Voltage
- Fully Portable
- Oil-free lubrication pump
- Square negative pressure dial
- Built-in bottle cavity
- Low noise
- Small size (H22CM x L29CM x D20CM)
- Easy to carry
- Fully portable unit incorporating three modes of electricity: 240V AC, 12V DC and Battery
- Suitable for vehicles, clinic, home and outdoor places

Technical Information:

Voltage:	AC 220V/240V 50HZ 12VDC
Power:	90 VA
Max Negative Pressure Value:	≥0.075 MPa
Pumping Rate:	≥15 Litres/Minute
Reservoir Volume:	1,000 mL
Noise:	≤65 Db (A)

CODE
SPBB

DESCRIPTION
LIBERTY 7E-D SUCTION PUMP 240V/12V/BATTERY 15L

AVAILABLE FROM:



MAJAC MEDICAL PRODUCTS Pty Ltd

Tel : 1300 138 578
Fax: 1300 138 612

ABN 20 086 942 421
www.majacmedical.com.au
Email: sales@majacmedical.com.au