

Dino-Lite Medical Digital Microscopes

Dermascope - Digital



Dino-Lite Medical Dermatoscopes are ideal for recording detailed examinations of the skin with accurate measurement & calibration functions. All models are TGA Listed and carry the CE Mark



FEDMEDM700P (DERMATOSCOPE) (Metal)

Resolution: 5 Mpixels
Magnification: 10x to 50x, 200x
Illumination: Built-in 8 white LED lights.
Image: Photograph, video or time lapsed
Dimensions: 10cm (Length), 3.2 cm (Diameter)
Contact Cap for Polarization: **Code P2**



AM413ZT (DERMATOSCOPE)

Resolution: 1.3 Mpixels
Magnification: 10x to 50x, 200x
Illumination: Built-in 8 white LED lights.
Image: Photograph, video or time lapsed
Dimensions: 10cm (Length), 3.2 cm (Diameter)
Contact Cap for Polarization: **Code P2**



The DinoLite Medical Digital Dermatoscope includes DinoCapture 2.0 software allowing the practitioner to capture pictures, record video, and to annotate images. This function includes measurement and calibration capabilities.

TEL: 1300 138 578 FAX: 1300 138 612 E-Mail: sales@majacmedical.com.au