

Dino-Lite Medical Digital Microscopes

Otoscope - Digital



OTOSCOPES

Digital handheld otoscopes are light weight and easy to use. High quality 1.3 Megapixel sensor provides the user the ability to view and record images of the ear canal/eardrums without loss of detail.



AMH-EUT (OTOSCOPE)

Resolution: 1.3 Mpixels
Magnification: 55x to 90x
Illumination: Built-in 6 white LED lights.
Image: Photograph, video or time lapsed
Dimensions: 11.9cm (Length), 3.2 cm (Diameter)
Weight: 95g
Specula: Disposable, SP3, SP4 and SP5



AM311H (OTOSCOPE)

Resolution: 640
Magnification: 20x to 30x
Illumination: Built-in 4 white LED lights.
Image: Photograph, video or time lapsed
Dimensions: 12cm (Length), 3.2 cm (Diameter)
Weight: 95g
Specula: Disposable, SP3, SP4 and SP5



Disposable Specula
3mm, 4mm and 5mm



TEL: 1300 138 578 FAX: 1300 138 612 E-Mail: sales@majacmedical.com.au