



DRESSINGS

The most comprehensive range of dressings to satisfy hospital, surgical & first aid requirements.

Low linting.
Highly adsorbent.
Available in sterile & non-sterile.
Comes in various sizes & lengths.



COMBINE PAD DRESSING (sterile)

Combine Dressings consists of a highly absorbent thick layer of fleece enclosed in a soft and comfortable non-woven fabric, pre-cut into manageable pieces. The non-woven cover pads have soft end seals to prevent pressure areas from forming. These are suitable for use, either in direct contact for medium to heavily exuding wounds, or as a secondary dressing.

RE-ORDER	DESCRIPTION	UNIT	CARTON QTY
NWC001	10cm x 10cm	50	16 packs x50
NWC002	10cm x 20cm	50	8 packs x50
NWC003	20cm x 20cm	50	5 packs x50



COMBINE ABSORBENT ROLL (non-sterile)

Combine Absorbent Rolls consist of a highly absorbent thick layer of fleece enclosed in a soft and comfortable non-woven fabric, folded over to form a flat tube. This roll can be cut into sizes that suit. These are suitable for use, either in direct contact for medium to heavily exuding wounds, or as a secondary dressing.

RE-ORDER	DESCRIPTION	UNIT	CARTON QTY
NWC100	10cm x 10m roll	10 rolls	10 rolls
NWC101	20cm x 10m roll	10 rolls	10 rolls



UNIVERSAL MULTI TRAUMA COMBINE DRESSINGS (sterile)

Universal Multi Trauma Combine Dressings consists of a highly absorbent thick layer of fleece enclosed in a soft and comfortable non-woven fabric, folded over to form a flat tube. The non-woven cover pads have soft end seals to prevent pressure areas from forming. These are suitable for use, either in direct contact for medium to heavily exuding wounds, as a secondary dressing or as a padding or protection for major wounds.

RE-ORDER	DESCRIPTION	UNIT	CARTON QTY
UD001	20cm x 91cm large	100	100
UD002	20cm x 43cm small	100	100



MAJAC MEDICAL PRODUCTS

"Let's Get It Right"

Sales Inquiries Tel:1300 138 578 Fax:1300 138 612
Business Office Tel: +61 7 3265 6355 Fax +61 7 3265 2729
e-Mail: sales@majacmedical.com.au
www.majacmedical.com.au

Your Complete Range of Cleaning & Infection Control Products

DRESSINGS



MEDIAM210 NON-ADHERENT DRESSINGS

Mediam210 consists of a highly absorbent pad which is heat bonded on both sides to a very thin highly embossed polyethylene film. Suitable for draining of exudate and reduction in trauma to the healing process.

RE-ORDER	DESCRIPTION	UNIT	CARTON QTY
NAD001	Mediam210 - 5cm x 5cm	50	10 boxes x50
NAD002	Mediam210 - 7.5cm x 7.5cm	50	10 boxes x50
NAD003	Mediam210 - 10cm x 10cm	50	8 boxes x50
NAD004	Mediam210 - 10cm x 20cm	50	8 boxes x50



BURNS DRESSINGS

Burn Dressings have a highly absorbent inner, laminated on both sides to a thin low adherent layer. Applied to wound sites to absorb exudate.

RE-ORDER	DESCRIPTION	UNIT	CARTON QTY
BD001	Burns dressing 75cm x 220cm	50	50
BD002	Burns dressing 75cm x 75cm	100	100



AsGUARD ISLAND DRESSINGS (sterile)

AsGuard Dressings combine an absorbent pad (Mediam210) with a soft and comfortable non woven backing layer for the simple and effective management of wound sites.

RE-ORDER	DESCRIPTION	UNIT	CARTON QTY
ISD001	AsGUARD - 7.5cm x 5cm	50	10 boxes x50
ISD002	AsGUARD - 8cm x 10cm	50	10 boxes x50
ISD003	AsGUARD - 15cm x 10cm	20	10 boxes x20
ISD004	AsGUARD - 20cm x 10cm	20	10 boxes x20



WOUND DRESSINGS (sterile)

These field wound dressings are a thick highly absorbent cotton pad attached to a gauze bandage. Ideal for first aid and emergency response.

RE-ORDER	DESCRIPTION	UNIT	CARTON QTY
WD001	No.13 small	12	48 dozen
WD002	No.14 medium	12	22 dozen
WD003	No.15 large	12	18 dozen



EYE PADS

Made from absorbent cotton and covered with non-woven rayon, these 6cm x 8cm eye pads are available in either individual sterile pre-packs or as non-sterile.

RE-ORDER	DESCRIPTION	UNIT	CARTON QTY
EP001	Eye pad 6cm x 8cm	50	24 boxes x50



MAJAC MEDICAL PRODUCTS

"Let's Get It Right"

Sales Inquiries Tel:1300 138 578 Fax:1300 138 612
Business Office Tel: +61 7 3265 6355 Fax +61 7 3265 2729
e-Mail: sales@majacmedical.com.au
www.majacmedical.com.au

Your Complete Range of Cleaning & Infection Control Products

AsGUARD CLEAR



AsGUARD Clear is an adhesive film dressing which has a high moisture vapour transfer rate. It provides an environment for moist wound healing and can also be used as a secondary dressing. This type of film dressing can be left on the wound site for extended periods thereby minimizing interference with the wound site. It also allows the natural healing process to occur while the excessive exudate evaporates.



RE-ORDER	DESCRIPTION	UNIT	CARTON QTY
ISB001	AsGuard Clear - 4cm x 5cm	50	10 boxes x50
ISB002	AsGuard Clear - 6cm x 7cm	50	10 boxes x50
ISB003	AsGuard Clear - 10cm x 12cm	25	10 boxes x25
ISB004	AsGuard Clear - 15cm x 20cm	25	10 boxes x25

AsGUARD PLUS+



AsGUARD PLUS+ combines a high moisture vapour transmission rate film and the cushioning and absorbency of a low adherent pad (Mediam210) together with an easy application system. The waterproof film allows the patients to shower or bath with it on without affecting the wound site or changing dressings. With the clear film backing the level of exudate can be easily monitored.



RE-ORDER	DESCRIPTION	UNIT	CARTON QTY
ISC001	AsGuard Plus - 4cm x 5cm	50	10 boxes x50
ISC002	AsGuard Plus - 6cm x 7cm	50	10 boxes x50
ISC003	AsGuard Plus - 10cm x 12cm	25	10 boxes x25
ISC004	AsGuard Plus - 15cm x 20cm	25	10 boxes x25

DISTRIBUTOR



MAJAC MEDICAL PRODUCTS

"Let's Get It Right"

Sales Inquiries Tel:1300 138 578 Fax:1300 138 612
 Business Office Tel: +61 7 3265 6355 Fax +61 7 3265 2729
 e-Mail: sales@majacmedical.com.au
www.majacmedical.com.au

Your Complete Range of Cleaning & Infection Control Products