

Sonica Ultrasonic Cleaners

Product

Description

Packing



Power and Technology in Ultrasonic Cleaning

SOLTEC designs and produces ultrasonic cleaners for professional use in all fields in which high cleaning standards must be achieved within the shortest possible timescale such as the electronic, pharmaceutical, dental, medical, military and aerospace sectors. Using ultrasonic energy and a process known as cavitation, ultrasonic cleaning is the most effective way of deep cleaning parts that are difficult to reach or complex different shaped items (mechanical parts and precision instruments).

In all these uses and applications, ultrasonic precision cleaning is the fastest and most economical way of effectively removing oil, grease, biofilm, organic matter, etc.



3 litres

Code: SOLSN2200M



4.5 litres

Code: SOLSN2400M

More than 30 models 3 to 30 litres

Code	Description	Tank Dimension
SOLSN2200M	2200M 3 ltr tank with basket and lid	240 x 140 x 100 mm
SOLSN2400M	2400M 4.5 ltr tank with basket and lid	300 x 150 x 100 mm
SOLSN3200M	3200M 6 ltr tank with basket and lid	300 x 240 x 100 mm
SOLSN3200LM	3200LM 6 ltr long tank with basket and lid	500 x 140 x 100 mm
SOLSN3300M	3300M 9.5 ltr tank with basket and lid	300 x 240 x 150 mm
SOLSN4200M	4200M 14 ltr tank with basket and lid	330 x 300 x 150 mm
SOLSN4300M	4300M 18 ltr tank with basket and lid	330 x 300 x 200 mm

Control Panel

Version M

The simplest version in the SONICA range is equipped with a timer which can be set for up to 15 minutes. It can also be set to manual operation

**Sonica has a range of Cleaners to include heaters, electronic controls and the new sweep system generator
For further information please contact our office.**

AVAILABLE FROM:



MAJAC MEDICAL PRODUCTS Pty Ltd

P.O. Box 74 Virginia QLD 4014

4/4 Johnstone Road, Brendale QLD 4500

Tel : 1300 138 578

Fax: 1300 138 612

www.majacmedical.com.au

Email: sales@majacmedical.com.au