

CliniVac®

Dental Aspiration & Suction Unit Detergent

CliniVac® is a Clinical Detergent Specially Formulated for use in Dental Aspiration Lines & Separation Systems

A Powerful Effective Suction Line Cleaner

- ✓ **Rapid and effective hard surface detergency**
Promotes penetration of the adhering bio-soil.
- ✓ **Promotes “wetting” of hydrophobic surfaces**
Ensures that the cleaning solution has contact with the entire lumen surface.
- ✓ **Corrosion inhibitor** incorporated into the product to protect metal surfaces.
- ✓ **Highly buffered** to effectively ionize heavy protein and bio-soil loads. This ensures that the washing solution maintains the optimal pH throughout the entire length of the tube, helping to prevent bio-film build-up.



CliniVac® is a product of Quality Medical Innovations Pty. Ltd.
CliniVac® is registered Trade Mark used under Licence.



MAJAC MEDICAL PRODUCTS

TELEPHONE: 1300 138 578 FACSIMILE: 1300 138 612

E-Mail: sales@majacmedical.com.au, www.majacmedical.com.au



Silver Partner