

Instrument Cleaning Brushes

Product

Description




These Cleaning Brushes are approx 9 cm length with grid type handle. The bristles are strong and hard wearing.

Stainless steel bristles are available for those with a special need.

Both models can be autoclaved.

Code: 3181 – Stainless Steel (3 per pack)

Code: 3182 – Nylon (3 per pack)



There is a wide range of cleaning brushes available to you. Tube brushed are available for cleaning lumens from ultra small to large that is the right size to clean cannulated instruments. Choose from extra long to small diameter or from the standard size tube brushed. All brushed are made of nylon bristles imbedded in twisted stainless steel wire.



These scrub brushes are durable and the fine bristles assist when cleaning contaminated instruments.

- Autoclavable scrub brush
- Black nylon fibres
- Surgeons choice

Code: BN12 (Box of 12)

AVAILABLE FROM:



MAJAC MEDICAL PRODUCTS Pty Ltd

Tel : 1300 138 578

Fax: 1300 138 612

ABN 20 086 942 421

www.majacmedical.com.au

Email: majac@optusnet.com.au