



# Fetal Doppler / Vascular Doppler

## Bestman Pocket Fetal Doppler



Built in Rechargeable Battery (lasts up to 6 hours continuous use)  
 Can use standard 9V battery  
 2Mhz probe standard (optional 5.0Mhz probe for umbilical blood stream)  
 Size 20.0cm x 16.5cm x 6.5cm (Only 270g )  
**Code: BESBF-500D**

## Bestman BF-560 Multi-purpose Doppler



LCD Bar Graph displays Flow Velocity  
 Head Phone Jack  
 Back Light On/Off  
 Large Clear Speaker with Adjustable Volume  
 Size: 35 x 90 x 130mm (only 500g)  
 Optional Vascular 8Mhz or Vascular 5Mhz Probe

### Standard Configuration Options

BESBF-560-F  
 BESBF-560-V

Foetal Doppler (2Mhz Probe)  
 Vascular Doppler (8Mhz Probe)



### Additional Probe Options

BESBF-560-2Mhz  
 BESBF-560-5Mhz  
 BESBF-560-8Mhz

Foetal Probe  
 Vascular Probe  
 Vascular Probe



## Bestman Vascular Doppler with Printer and LCD Display

8Mhz or 5Mhz Vascular Probe  
 LCD Screen displays Flow Velocity Peak, Mean, Heart Rate and Time  
 Adjustable Volume and Time/Date Display  
 Built in High Precision Printer with 55mm paper  
 RS-232 Serial output for Computer Connection  
 Weight 1.4Kg  
**Code: BESBV-620VP**



# MAJAC MEDICAL PRODUCTS

*"Let's Get It Right"*

Sales Inquiries Tel:1300 138 578 Fax:1300 138 612  
 Business Office Tel: +61 7 3265 6355 Fax +61 7 3265 2729  
 e-Mail: sales@majacmedical.com.au  
 www.majacmedical.com.au

Your Complete Range of Cleaning & Infection Control Products