

Standard Exam Couch PLUS with braked castors and gas strut backrest

Code: 383.161.SP \$2,775.00



The Standard Exam Couch PLUS (+ electric height and Gas strut back rest) is a versatile examination couch.

Specifications:

- * Grey base / Balmain Navy upholstery / 710mm wide
- * Overall length 1930mm / backrest section 760mm [+80°]
- * Electric height adjustment [45 - 95cm] with foot control
- * Two leg / swing arm direct lifting system with bearings on moving parts.
- * Single action retractable castors.
- * Lifting capacity 200 kg

Condor brand medical grade heavy-duty vinyl fabric, delivers excellent strength and wear characteristics and is fire retardant, perspiration and oils resistant. A choice of colours, refer brochure). Note: Electric back rests are available and extra. TGA approved ARTG No.123573.

Treatment Chair / Couch D100 - 3 section

Code: 720.581.0 standard \$3,345.00

Code: 720.681.0 with electric leg rest \$4,340.00



Standard Specifications

- * Grey base / TBA upholstery / 600mm wide / Swing back arm rests
- * Four sections: 890mm contoured [+80°] / 465mm [reclines 10° when backrest lifted] / 310mm [+20/-50°] / 100mm leg extension plus two arm rests
- * Electric height and integrated backrest / seat tilt adjustment with foot control.
- * Single leg lifting system with bearings on moving parts
- * Single action retractable castors.
- * Lifting Capacity 175kg

Standard Accessories

Retractable castors for on-the-ground stability in use.
Swing back arm rests
Extending leg rest

Treatment Arm Board

Code: MAATB \$200.00



Treatment Arm Board Option - Accessory for Examination Couches.
Suitable for most Examination Couches.

The Treatment arm board comprises an upholstered arm-rest (400mm long, tapering from 150mm to 100mm wide) which can be located on either side of the couch.

It can rotate 360°, is height adjustable to 300mm and can be positioned up to 150mm from the edge of the couch. Factory fitted usually at time of purchase of couch.

Gynae Knee Supports (lithotomy)

Code: MAGKS \$475.00



Lithotomy Knee Supports - Standard on Gynae Exam Couches - A Factory fitted option on other Exam Couches.

The black plastic supports are specially moulded for excellent support and comfort. They can be positioned in a wide range of angles to suit the particular patient and procedure using the universal joint.

They can rotate 360°, are height adjustable to 200mm and can be positioned up to 150mm from the edge of the couch. The supports can be removed when not required. They are a factory fitted accessory, usually prior to delivery of couch.

Side Rails

Code: MADSR Cot Side Rail (71cm) \$125.00 each

Code: MADSRP Upholstered Side Rail (95x40 cm) \$220.00 each



For use, the cot side is lifted up and drops into the support brackets. When not in use the cot side folds down.

Left and right side available.

Gynaecological Exam Couch

Exam G Gynaecological Examination Couch

Code: 6100.161.0 \$3,415.00



The Exam G Gynaecological Examination Couch (Table) - height adjustable - Gynaecological Couch. (electric height and gas strut back rest). The Gynae is a versatile couch suitable for either general or gynaecological examination and has a perineal recess cut away in the foot section. **The foot supports are inclusive and under knee (lithotomy) supports available as an optional extra.**

Specifications:

- * Grey base / Balmain Navy upholstery / 710mm wide
- * Overall length 1930mm / backrest section 760mm [+80°]
- * Electric height adjustment [45 - 95cm] with foot control
- * Two leg / swing arm direct lifting system with bearings on moving parts.
- * Single action retractable castors.
- * Lifting capacity 200kg

Condor brand medical grade heavy-duty vinyl fabric, delivers excellent strength and wear characteristics.

A choice of colours, (refer brochure). Note: Electric back rest is available and extra. TGA approved ARTG No.123573.

Exam G Gynaecological Examination Couch, Twin Back Rest

Code: 611.161.GT \$3,770.00



The Exam G Gynaecological Examination Couch (Table) - height adjustable - "GT" TWIN BACK REST Gynaecological Couch. The "GT" has a second backrest to provide support during gynaecological examinations and a perineal recess cut away in the foot section. (electric height and Gas strut back rests) This is a versatile couch suitable for either general or gynaecological examination. **Foot supports are standard and Knee (lithotomy) supports are an optional accessory.**

Specifications:

- * Grey base / Balmain Navy upholstery / 710mm wide
- * Overall length 1930mm / backrest section 760mm [+80°]
- * Electric height adjustment [45 - 95cm] with foot control
- * Two leg / swing arm direct lifting system with bearings on moving parts.
- * Single action retractable castors.
- * Lifting capacity 200kg

Condor brand medical grade heavy-duty vinyl fabric, delivers excellent strength and wear characteristics.

A choice of colours, (refer brochure). Note: Electric back rest is available and extra. TGA approved ARTG No.123573.

Gynae Foot Support Pads

Code: MAGFP \$360.00



Gynae Foot Support Pads are an accessory for the range of Examination Couches. The pads are adjustable, easy to use and comfortable for the patient. The Supports can be removed when not required. They are a factory fitted accessory, usually prior to delivery of the couch.

Upholstery Colour Options:

(Colours are approximate only)



Condor brand medical grade heavy-duty vinyl fabric. The exclusive formulation delivers a superior softness for comfort including excellent strength and wear characteristics. Formulated to inhibit hardening and breakdown of the PVC surface due to perspiration and oils; light fast; fire retardant; sanitised; stain resistant.

Maintenance

Recommended cleaning with chemically neutral soap and warm water, then wipe clean. Do not dry-clean or use abrasive or solvent based products.

Gas Lift Stools



Standard Specifications

* 40cm diameter stool that adjusts between approx. 40 and 55cm by gas strut.
Upholstery colour to match your couch or chair.

27550.XX0

Gas lift stool with back \$285.00

27500.XX0

Gas lift stool without back \$205.00

Equipment Trolleys



ET602/ET605 Trolley \$160.00



ET602/ET605 Trolley with drawers \$205.00

Ideal Trolley for ECGs, Spirometers, Resuscitators & Defibrillators

Standard Specifications

- * Sturdy black tubular steel trolleys with grey steel shelves
- * Overall dimensions: ET602 - 45cm w x 39cm d x 72cm height
ET605 - 45cm w x 39cm d x 90cm height

Black frame with two grey steel shelves 40 x 39cm. Supplied flat for DIY assembly. Optional drawers. White Drawers. Internal dimensions: 350mm wide x 325mm deep x 110mm high. One of the shelves can be mounted in the opposite direction to the other to give improved access and visibility to the lower shelf.

71201 ET602 Equipment Trolley (72 cm height)

71203 ET605 Equipment Trolley (90 cm height)

71202 ET602 Equipment Trolley Drawer (72 cm height)

71204 ET605 Equipment Trolley Drawer (90 cm height)



MAJAC MEDICAL PRODUCTS

TELEPHONE: 1300 138 578 FACSIMILE: 1300 138 612

E-Mail: sales@majacmedical.com.au, www.majacmedical.com.au



Silver Partner