

Examination Couch (Table) - Height Adjustable

Code: 853.181.2

“GP” Examination Couch (Table) – Height Adjustable, Designed for General Practice

Meets the minimum specifications of disability advocates (www.racgp.org.au/standards/513)

Code: 853.181.2 with **manual backrest** \$1,790.00

Code: 853.181.G2 with **gas strut backrest** \$1,895.00



The “GP” Exam Couch - (electric height and manual backrest). This model is a versatile examination couch.

Specifications:

- * Grey base / Balmain Navy upholstery / 710mm wide
- * Overall length 1930mm / backrest section 680mm [0/+80°]
- * Electric height adjustment [45 - 95cm] with foot control
- * Two leg / swing arm direct lifting system with bearings on moving parts.
- * Single action retractable castors.
- * Lifting capacity 200 kg

Grey Base / Balmain Navy upholstery. Condor brand medical grade heavy-duty vinyl fabric, delivers excellent strength and wear characteristics.

Formulated to inhibit hardening and breakdown of the surface due to perspiration and oils; fire retardant; sanitised; stain resistant.

Manual back rest. Breathing Hole and plug is a factory fitted option for additional \$75.00. TGA approved ARTG No.123573.

Note: Breathing Hole extra

code: 851.181.2

“GP” Examination Couch (Table) – Height Adjustable, Designed for General Practice 3 Section

Meets the minimum specifications of disability advocates (www.racgp.org.au/standards/513)

\$1950.00



Examination Couch – Height Adjustable – 3 Section.

Specifications:

- * Grey base / Balmain Navy upholstery / 710mm wide
- * Three sections: 360mm with breathe hole & plug (+50 / -80°) / 560mm / 940mm (+80deg) – Overall length 1930 mm
- * Electric height adjustment (50 - 100cm) with foot control
- * Two leg / swing arm direct lifting system with bearings on moving parts.
- * Single action retractable castors.
- * Lifting capacity 200kg
- * One breathe hole & plug is standard

Code: MAATB

Treatment Arm Board

\$175.00



Treatment Arm Board Option - Accessory for Examination Couches.

Suitable for most Examination Couches.

The Treatment arm board comprises an upholstered arm-rest (400mm long, tapering from 150mm to 100mm wide) which can be located on either side of the couch.

It can rotate 360°, is height adjustable to 300mm and can be positioned up to 150mm from the edge of the couch. Factory fitted usually at time of purchase of couch.

Code: MAGKS

Gynae Knee Supports (lithotomy)

\$475.00



Lithotomy Knee Supports - Standard on Gynae Exam Couches - A Factory fitted option on other Exam Couches.

The black plastic supports are specially moulded for excellent support and comfort. They can be positioned in a wide range of angles to suit the particular patient and procedure using the universal joint.

They can rotate 360°, are height adjustable to 200mm and can be positioned up to 150mm from the edge of the couch. The supports can be removed when not required. They are a factory fitted accessory, usually prior to delivery of couch.

Code: MADSR & MADSRP

Side Rails

Code: MADSR \$125.00 each

Padded Side Rails, Code: MADSRP \$220.00 each



Factory fitted option for Exam Couches

For use, the cot side is lifted up and drops into the support brackets. When not in use the cot side folds down.

Left and right side available.

Code: 383.161.SP

Standard Exam Couch PLUS with braked castors and gas strut backrest

\$2,775.00



The Standard Exam Couch PLUS (+ electric height and Gas strut back rest) is a versatile examination couch.

Specifications:

- * Grey base / Balmain Navy upholstery / 710mm wide
- * Overall length 1930mm / backrest section 760mm [+80°]
- * Electric height adjustment [45 - 95cm] with foot control
- * Two leg / swing arm direct lifting system with bearings on moving parts.
- * Single action retractable castors.
- * Lifting capacity 200 kg

Condor brand medical grade heavy-duty vinyl fabric, delivers excellent strength and wear characteristics and is fire retardant, perspiration and oils resistant.

A choice of colours, refer brochure). Note: Electric back rests are available and extra. TGA approved ARTG No.123573.



MAJAC MEDICAL PRODUCTS

TELEPHONE: 1300 138 578 FACSIMILE: 1300 138 612

E-Mail: sales@majacmedical.com.au, www.majacmedical.com.au



Silver Partner

Code: 6100.161.0

Exam G Gynaecological Examination Couch

\$3,415.00



The Exam G Gynaecological Examination Couch (Table) - height adjustable - Gynaecological Couch. (electric height and gas strut back rest). The Gynae is a versatile couch suitable for either general or gynaecological examination and has a perineal recess cut away in the foot section. The foot supports are inclusive and under knee (lithotomy) supports available as an optional extra.

Specifications:

- * Grey base / Balmain Navy upholstery / 710mm wide
- * Overall length 1930mm / backrest section 760mm [+80°]
- * Electric height adjustment [45 - 95cm] with foot control
- * Two leg / swing arm direct lifting system with bearings on moving parts.
- * Single action retractable castors.
- * Lifting capacity 200kg

Condor brand medical grade heavy-duty vinyl fabric, delivers excellent strength and wear characteristics.

A choice of colours, (refer brochure). Note: Electric back rest is available and extra. TGA approved ARTG No.123573.

Code: 611.161.GT

Exam G Gynaecological Examination Couch, Twin Back Rest

\$3,770.00



The Exam G Gynaecological Examination Couch (Table) - height adjustable - "GT" TWIN BACK REST Gynaecological Couch. The "GT" has a second backrest to provide support during gynaecological examinations and a perineal recess cut away in the foot section. (electric height and Gas strut back rests) This is a versatile couch suitable for either general or gynaecological examination. Foot supports are standard and Knee (lithotomy) supports are an optional accessory.

Specifications:

- * Grey base / Balmain Navy upholstery / 710mm wide
- * Overall length 1930mm / backrest section 760mm [+80°]
- * Electric height adjustment [45 - 95cm] with foot control
- * Two leg / swing arm direct lifting system with bearings on moving parts.
- * Single action retractable castors.
- * Lifting capacity 200kg

Condor brand medical grade heavy-duty vinyl fabric, delivers excellent strength and wear characteristics.

A choice of colours, refer brochure). Note: Electric back rest is available and extra. TGA approved ARTG No.123573.



MAJAC MEDICAL PRODUCTS

TELEPHONE: 1300 138 578 FACSIMILE: 1300 138 612

E-Mail: sales@majacmedical.com.au, www.majacmedical.com.au



Silver Partner

Gynae Foot Support Pads**\$360.00**

Gynae Foot Support Pads are an accessory for the range of Examination Couches. The pads are adjustable, easy to use and comfortable for the patient. The Supports can be removed when not required. They are a factory fitted accessory, usually prior to delivery of the couch.

Upholstery Colour Options:

(Colours are approximate only)



Condor brand medical grade heavy-duty vinyl fabric. The exclusive formulation delivers a superior softness for comfort including excellent strength and wear characteristics.

Formulated to inhibit hardening and breakdown of the PVC surface due to perspiration and oils; light fast; fire retardant; sanitised; stain resistant.

Maintenance

Recommended cleaning with chemically neutral soap and warm water, then wipe clean.

Do not dry-clean or use abrasive or solvent based products.