

Nitrile Examination Gloves

Product

Description



Description

Whilst Healthcare wearers require protection from exposure to infectious substances and cross-contamination, latex sensitivity and allergic reactions can be a concern. Nitrile gloves offer protection without the risk of latex reactions.

Features & Benefits

- Superior puncture and tear resistance
- Excellent tactile sensitivity and dexterity
- Textured fingertip and palm to enhance wet and dry grip
- Lavender colour provides quick, visual differentiation
- Beaded cuff to prevent roll down and for easier donning
- Power-free
- Latex-free
- Length 24cm
- Meets and exceeds Australian Standard requirements AS/NZ 4011:1997
- Manufactured in accordance with Quality System ISO 9002

Indications

Where latex-free protection is required

Code	Description	Item per unit
36361340	Small	Box/100
36361341	Medium	Box/100
36361342	Large	Box/100
36361343	X-Large	Box/100

AVAILABLE FROM:



MAJAC MEDICAL PRODUCTS Pty Ltd

ABN 20 086 942 421

Tel : 1300 138 578

Fax: 1300 138 612

www.majacmedical.com.au
Email: sales@majacmedical.com.au

Gloves - Procedure

Product


Description



Nitrile and Latex Gloves

All Sizes Available: Large, Medium, Small, Extra Small

Non Sterile – Ambidextrous - Textured



Nitrile Gloves

Nitrile gloves are excellent as a barrier and protection for medical procedures. The softness and stretch-ability ensures ultimate superiority and comfort.

Physical Features

- Cuff: Beaded
- Texture: Full textured & finger textured
- Colours: White, dark blue and light blue

Material Compounds

Our nitrile gloves are free from the following chemicals:

- Plasticiser
- Nonyl Phenol
- Styrene
- Latex Protein

Code: ADGN

Unique Features

- Excellent alcohol resistance
- Outstanding stress relaxation
- Excellent durability
- Extensible comfort
- Good hydration properties
- Water durability
- Elasticity better than Vinyl
- Excellent puncture resistance
- Good resistance to chemicals
- Can be sterilised by Gamma Radiation

Latex Gloves

Both use powdered and powder-free gloves are formulated for use during medical procedures. Powdered latex gloves were produced to meet the recommended allowable protein level.



Material Compounds

Our natural rubber latex gloves are free from the following chemicals:

- Plasticiser
- Acrylonitrile
- Nonyl Phenol
- Styrene

Code: ADGLPF – Latex Power Free

Code: ADGLP – Latex Powdered

Unique Features

- Excellent alcohol resistance
- Outstanding stress relaxation
- Excellent gripping properties
- Extensible comfort
- Water durability
- Excellent elasticity
- Can be sterilised by Gamma Radiation

AVAILABLE FROM:



MAJAC MEDICAL PRODUCTS Pty Ltd

P.O. Box 74 Virginia QLD 4014

4/4 Johnstone Road, Brendale QLD 4500

Tel : 1300 138 578

Fax: 1300 138 612

www.majacmedical.com.au

Email: sales@majacmedical.com.au