

CENTURION MINOR SURGERY LIGHTING

The Centurion is a revolutionary mid-range surgery light which offers the performance and features found on higher-priced surgery lights. The Centurion includes a highly efficient, single reflector design, assuring optimal shadow reduction, and the 21st Century Series exclusive Centra arm delivery system.



- Prefocused spot at a distance of 36" (914mm) from lamphead providing homogenous light with 4,000 ft candles (43,000lux) and color temperature of approximately 4,000° Kelvin
- Easily replaceable, low cost, long life 12 volt/90 watt halogen bulb
- 15" (380mm) diameter elliptical reflector
- IR absorbing, color-correcting glass cylinder
- 28" (711mm) articulating arm with integral compression spring
- 23" Extension Arm
- Arm rotates 40° up and down from horizontal plane
- Lamp head and arm both rotate 540°
- Removable sterilizable handle
- Additional positioning handle mounted on perimeter
- Available in portable floor, wall, single or dual ciling mounts, combination ceiling mounts with Centura or Centry and ChuttleTrak ceiling track mount.



Wall model 021512



Single Ceiling model 021514



Portable Floor model 021513



MAJAC MEDICAL PRODUCTS

"Let's Get It Right"

Sales Inquiries Tel:1300 138 578 Fax:1300 138 612
Business Office Tel: +61 7 3265 6355 Fax +61 7 3265 2729
e-Mail: sales@majacmedical.com.au
www.majacmedical.com.au

Your Complete Range of Cleaning & Infection Control Products

CENTURION eXceL SURGERY LIGHT

Centurion eXceL represents the new standard in minor surgery lighting, featuring a larger reflector, higher output, a larger pattern and improved shadow control. Available in all of the normal 21st Century Series mounting options: Wall, Portable, Single, Dual or Combo Fixed Ceiling Mounts.



- Prefocused 6"(150mm) spot, at 1 meter distance
- 6,000 ft candles (64,500 lux) output at 1 meter
- 16"(406 mm) diameter, multi-faceted reflector for excellent shadow control
- IR absorbing, color-correcting glass cylinder
- 90 watt halogen bulb - quick bulb replacement
- 28" (711 mm) articulating arm with integral compression spring
- 23" extension arm
- Lamp head and arm both rotate 540 degrees
- Arm moves 40 degrees up and down from horizontal plane
- Standard removeable, sterilizable handle
- Additional positioning handle on perimeter of light head
- Removable, sterilizable handle
- Available Mountings: Wall, Portable, Single, Dual or Combo fixed ceiling mounts, Chuttletrak ceiling track mounts



Wall model 051512



Single Ceiling model 051514



Portable Floor model 051513



MAJAC MEDICAL PRODUCTS

"Let's Get It Right"

Sales Inquiries Tel:1300 138 578 Fax:1300 138 612
Business Office Tel: +61 7 3265 6355 Fax +61 7 3265 2729
e-Mail: sales@majacmedical.com.au
www.majacmedical.com.au

Your Complete Range of Cleaning & Infection Control Products