



MAJAC MEDICAL PRODUCTS Pty Ltd

P.O. Box 74 Virginia QLD 4014
11/209 Robinson Road
Geebung QLD 4034
Australia
www.majacmedical.com.au

ABN 20 086 942 421

Tel: 1300 138 578


Fax: 1300 138 612

Email: sales@majacmedical.com.au
admin@majacmedical.com.au

MDScope Quick Installation Guide

Step 1. Hardware Installation

Please follow the pictorial guide below for input / output connections:



When you first connect the USB Capture to your computer a Hardware Assistant will appear. Select "**Cancel**" button to close "Hardware Update Wizard" window.



Step 2. Insert the "Installations CD" into your CD drive

The CD automatically displays the AutoRun windows if Autorun is enable in your computer. Select "**Install Software**" button to complete the Driver & Applicatron installation.

