

# Finger Tip Pulse Oximeter

## Product

## Description

Now Oximetry has become the 3rd most important medical practice or homecare instrument following Blood Pressure Meter and Glucose Meter.

Pulse oximeter is a very useful and efficient device to show patient's Haemoglobin Saturation condition, greatly prevents patient from threat of respiration disease, chest disease, etc.

With low cost, reliability and extreme light weight portability, fingertip pulse oximeter can be afforded by every family, like a blood pressure meter & glucose meter.

## MD300C21

### Fingertip Tip Pulse Oximeter - Adult

low power consumption and automatic power-off



This small MD300C21 oximeter monitors the patient's oxygen saturation and pulse-rate (PR) and visually displays the recorded SpO2 percentage and PR values on a bright LCD screen. For added convenience, the direction of the display can be rotated 180° for easy overview of data. This light and compact oximeter is easy to use and operates for approximately 30 hours with two AAA batteries. Indication of low power consumption and automatic "power off" eight seconds later. Suitable for the clinic, hospital or home use.



## MD300C5

### Pulse Oximeter Child Finger LCD



WORLD FIRST! Designed for Children 4 weeks to 12 years of age  
OLED display; Finger tip sensor  
Parameter: SpO2, PR, Pleth Bar – Auto Power Off  
2 AAA battery power supply;  
approximately 30 hours battery life



AVAILABLE FROM:



**MAJAC MEDICAL PRODUCTS Pty Ltd**

P.O. Box 74 Virginia QLD 4014

11/209 Robinson Road East, Geebung QLD 4034

Tel : 1300 138 578

Fax: 1300 138 612

[www.majacmedical.com.au](http://www.majacmedical.com.au)

Email: [sales@majacmedical.com.au](mailto:sales@majacmedical.com.au)