

# Finger Tip Pulse Oximeter

## Product

## Description

Now Oximetry has become the 3rd most important medical practice or homecare instrument following Blood Pressure Meter and Glucose Meter.

Pulse oximeter is a very useful and efficient device to show patient's Haemoglobin Saturation condition, greatly prevents patient from threat of respiration disease, chest disease, etc.

With low cost, reliability and extreme light weight portability, fingertip pulse oximeter can be afforded by every family, like a blood pressure meter & glucose meter.

## MD300D

### Fingertip Tip Pulse Oximeter

low power consumption and automatic power-off



This small MD300D oximeter monitors the patient's oxygen saturation and pulse-rate (PR) and visually displays the recorded SpO2 percentage and PR values on a bright LCD screen. For added convenience, the direction of the display can be rotated 180° for easy overview of data. This light and compact oximeter is easy to use and operates for approximately 30 hours with two AAA batteries. Indication of low power consumption and automatic "power off" eight seconds later. Suitable for the clinic, hospital or home use.



Price \$285.00 (ex GST)

<p><b>Patient Range</b> Adult, Pediatric</p> <p><b>SpO2</b> Range 0% - 99% Resolution 1% Accuracy 70%-100%: +/- 2% 0%-69%: unspecified</p> <p>Method Dual wavelength LED</p> <p><b>Heart (Pulse)Rate</b> Range 30-240 bpm Resolution 1 bpm Accuracy +/- 2 bpm or +/-2%</p>	<p><b>Display</b> Type LCD display 7 display modes</p> <p>Parameters SpO2, PR, Pleth bar. Waveform Low battery</p> <p><b>Information</b> Environmental 5°C - 40°C Operating Temperature -10°C - 40°C Storage Temp Operating Humidity 15% RH – 95% RH Storage Humidity 10% RH – 95% RH</p> <p><b>Classification</b> Type of protection internally powered equipment Safety IEC Standard 61601-1-1</p>	<p><b>Mechanical</b> Size 58 x 32 x 34 mm</p> <p>Weight 50g (including 2 batteries)</p> <p><b>Battery</b> Type 2 AAA size 1.5V Alkali batteries or 2 AAA 1.5V rechargeable batteries</p> <p>Operation time minimum 30 hours for normal operation</p>
--	--	--



AVAILABLE FROM:



**MAJAC MEDICAL PRODUCTS Pty Ltd**

P.O. Box 74 Virginia QLD 4014

11/209 Robinson Road East, Geebung QLD 4034

Tel : 1300 138 578

Fax: 1300 138 612

[www.majacmedical.com.au](http://www.majacmedical.com.au)

Email: [sales@majacmedical.com.au](mailto:sales@majacmedical.com.au)