

Handheld Pulse Oximeter

Product

Description

MD300A

Handheld Pulse Oximeter

MMED6000RD-MD300A is a specially designed Handheld Pulse Oximeter with SpO2 wave display designed for neonatal, pediatric and adult patients. Its compact design enables MD300A a wide use within clinical, homecare and emergency care environments

Pulse Oximeter HandHeld with Large Colour LCD
Large color LCD display; Handheld
Parameter: SpO2, PR, Pleth waveform, Pleth Bar
Audio and Visual alarm

DC or battery power supply

24 hours data storage & browser, SpO2 & PR trend.

Suitable for adult, pediatric and neonatal patient

Long battery life: approximately 20 hours

AC/DC adapter optional



Now available child (Pediatric) size sensor, for model MD300A and MD300I



MD Paed Y Sensor

Child Reusable Sensor Y Type
Suitable for Child 12 weeks old+
includes adhesive tape

For use with MD300A and MD2000A Pulse Oximeters

MD Child Finger Probe

Child Reusable Sensor (same as the Adult sensor but smaller)
Suitable for Child 12 weeks old+

For use with MD300A and MD2000A Pulse Oximeters



AVAILABLE FROM:



MAJAC MEDICAL PRODUCTS Pty Ltd

P.O. Box 74 Virginia QLD 4014

11/209 Robinson Road East, Geebung QLD 4034

Tel : 1300 138 578

Fax: 1300 138 612

www.majacmedical.com.au

Email: sales@majacmedical.com.au