

Handheld Pulse Oximeter

Product

Description

MD300I

Handheld Pulse Oximeter

Special Features

- Displays SpO2, PR, SpO2 waveform (pleth or line)
- PR tone modulation function
- Displays signal strength
- 10 patients ID setup & 72 hours recording storage
- Transmit data to PC (with optional software)
- Three display modes
- Alarm information: Probe off, Finger out, Low power
- Low power consumption and up to 24 hours continuous operation
- Light & compact design makes it convenience to carry
- Automatic power-off
- Low voltage indicator, Low power consumption, 2 AA alkaline or rechargeable batteries



Price \$649.00 (ex GST)

Technical Specifications

SpO2 Range 0% - 100% Resolution 1% Accuracy 80%-99%: +/- 2% 70%-80%: +/- 3% 0% - 69%: unspecified	Record Patient ID 10 patients Date record Up to 72 hours	Safety IEC Standard 61601-1-1
Heart (Pulse) Rate Range 30-254 bpm Resolution 1 bpm Accuracy +/- 2 bpm or +/- 2%	Data transmission Transmission method Data cable Interface Oximeter side: DB9 PC side: USB	Mechanical Size 130 x 40 x 27 mm Weight 120g (including batteries)
Alarm Alarm: Probe off, Finger out, Low power Modes Visual information	Environmental Operating Temperature 0°C - 50°C Storage Temp -10°C - 60°C Operating Humidity 15% RH - 95% RH Storage Humidity 10% RH - 95% RH	Power Type 2 AA Alkaline or Rechargeable batteries Operation time about 24 hours for normal operation
Display Type LCD Parameters SpO2, PR, Pleth Waveform, Pleth bar Mode: 3 display modes	Classification Type of protection internally powered equipment	



Optional Data Output Software for MD300I and MD300W Price: \$199.00 (ex GST)

AVAILABLE FROM:



MAJAC MEDICAL PRODUCTS Pty Ltd

P.O. Box 74 Virginia QLD 4014
 11/209 Robinson Road East, Geebung QLD 4034

Tel : 1300 138 578

Fax: 1300 138 612

www.majacmedical.com.au

Email: sales@majacmedical.com.au