

Wrist Mounted Oximeter

Product

Description

MD300W

Wrist Mounted Oximeter

Special Features

- Sleep monitoring
- Wrist measurement, convenient clinic operation
- Compact and lightweight design
- Application from adult to neonatal
- Low power consumption, working with 1 AAA alkaline battery
- Data transmission, storage and replay
- 1-10 patient ID setting
- Low voltage indicator



Price \$949.00 (ex GST)

Sleep Monitoring



Monitor during sleeping



Wrist measurement



Data transmission

<p>Patient Range Adult, Pediatric and Neonatal</p> <p>SpO2 Range 0% - 99% Resolution 1% Accuracy 80%-99%: +/- 2% 70%-80%: +/- 3% 0%-69%: unspecified Dual</p> <p>Heart (Pulse) Rate Range 0-254 bpm Resolution 1 bpm Accuracy +/- 2 bpm or +/- 2%</p> <p>Alarm Alarm Lower power indicator, probe off, finger off</p>	<p>Display Type LCD</p> <p>Parameters SpO2, PR, bargraph</p> <p>Environmental Operating Temperature 5°C - 40°C Storage Temp -10°C - 40°C Operating Humidity 15% RH - 80% Storage Humidity 10% RH - 80% RH</p> <p>Data transmission Transmission method data cable</p>	<p>Battery Type 1 x AAA Alkaline batteries Operation time about 12 hours for normal operation</p>
---	--	--

Optional Data Output Software for MD3001 and MD300W Price: \$199.00 (ex GST)

AVAILABLE FROM:



MAJAC MEDICAL PRODUCTS Pty Ltd

P.O. Box 74 Virginia QLD 4014

11/209 Robinson Road East, Geebung QLD 4034

Tel : 1300 138 578

Fax: 1300 138 612

www.majacmedical.com.au

Email: sales@majacmedical.com.au