

## Bowie-Dick Test Set

For Hospital Sterilisers or Sterilisers over 60 Litres  
(For both Hollow & Porous Loads)



### Daily Air Removal/Steam Penetration Test

Efficient removal of air is a critical parameter for the proper function of pre-vacuum (Class B) sterilisers. A steam penetration test should be performed in an empty steriliser at the beginning of every day to document the efficiency of the vacuum pumps and to indicate possible leaks and /or poor steam quality.

Accurate, sensitive and reproducible results warning of possible air leaks with 6 test pads to indicate pass or failure.

- ✓ Validation according to test pack DIN EN 285
- ✓ Best interpretation of colour change from **YELLOW** to **BLACK**. Non reversible indicator.
- ✓ Ink is free of any toxic components
- ✓ Colour Control Indicator has adhesive tape for attaching directly to sterilizer record system. eg The Trakka Sterilizer Record & Tracking System.

**Code: 15.100.072** *Bowie-Dick Test Set for Hollow & Porous loads in hospital pre-vacuum steam sterilisers. Includes Stainless Steel Helix device, 100 indicator strips, cotton bag & instructions for use.*

**Code: 15.100.050** *Bowie-Dick Test Refill Set for 15.100.072. Includes 250 indicator strips, 1 washer, 1 silicon junction, cotton bag & instructions for use.*



**MAJAC MEDICAL PRODUCTS**

TELEPHONE: 1300 138 578 FACSIMILE: 1300 138 612

E-Mail: [sales@majacmedical.com.au](mailto:sales@majacmedical.com.au), [www.majacmedical.com.au](http://www.majacmedical.com.au)



Silver Partner