

M+PACK – Self Sealing Sterilisation Pouch

Product

Description

M+PACK is a hospital grade self seal sterilisation pouch that conforms to the European standard EN 868-3.

The features include a pre-fold for accurate and secure closure this feature guarantees a complete air tight seal each time, dual process indicators, strong triple rail seals, thumb notch for ease of opening, batch numbered for traceability, plain paper backing for superior steam transmission with no ink contamination and a thumb notch to facilitate a clean peel when opening.

The fold and glue strip is 17mm wide and the perforated pre fold ensures an air tight seal. The protective strip over the fold may be removed and placed inside the pouch to act as an additional Process Indicator, if required.

M +PACK The perfect self seal Sterilising Pouch



Note: Air tight seal at fold (No Gaps) and protective strip “indicator” inside the pack and the steam process indicator on the outside of the pack (bottom left). Thumb notch for easy opening.

Re-Order	Description	Price
1135255	Pouch - Self Sealing Pouch - 133 x 255 (Box of 200)	25.00
1190330	Pouch - Self Sealing Pouch - 190 x 330 (Box of 200)	71.80
1230355	Pouch - Self Sealing Pouch - 230 x 350 (Box of 200)	91.80
1057100	Pouch - Self Sealing Pouch - 57 x 100 (Box of 200)	9.20
1070230	Pouch - Self Sealing Pouch - 70 x 230 (Box of 200)	13.30
1090135	Pouch - Self Sealing Pouch - 90 x 135 (Box of 200)	11.70
1090230	Pouch - Self Sealing Pouch - 90 x 230 (Box of 200)	14.80

AVAILABLE FROM:



MAJAC MEDICAL PRODUCTS Pty Ltd

ABN 20 086 942 421

Tel : 1300 138 578

Fax: 1300 138 612

www.majacmedical.com.au
Email: sales@majacmedical.com.au