

Resuscitator

Product

Description

MR100AD



Features

- Adult resuscitator 1500ml with 60cm H2O pop-off valve
- 2500ml Oxygen reservoir
- Two durable adult masks #5145, #5144
- Two guedel airways
- See through carry case with instruction manual

Code	Description
MR100AD	Durable Silicon Adult Resuscitation Set

MR100CH



Features

- Child resuscitator 550ml with 40cm H2O pop-off valve
- 2500ml Oxygen reservoir
- Two child masks 5143, 5122
- Two guedel airways
- See through carry case with instruction manual

Code	Description
MR100CH	Durable Silicon Child Resuscitation Set

MR100IN



Features

- Infant resuscitator 280ml with 40cm H2O pop-off valve
- 600ml Oxygen reservoir
- Two infant masks 5121, 5120
- Two guedel airways
- See through carry case with instruction manual

Code	Description
MR100IN	Durable Silicon Infant Resuscitation Set

1325



Features

- Bacterial air filter, Oxygen port
- Leakage free air cushion mask, silicon duck bill one way valve
- Detachable mouthpiece - mask can then be connected to other resuscitation devices
- **Specification:**
- Inspiratory Resistance : 4.5cm H2O
- Expiratory Resistance : 3.5cm H2O
- Material : PVC - Single patient use
- Dead Space : 125ml Connections: 15 mm OD
- Oxygen Concentration : 40-80% (depending on flow and ventilation cycles)
- Measurements (case) : 132 x 107 x 46mm

Code	Description
1325	Pocket Resuscitator with Filter

AVAILABLE FROM:



MAJAC MEDICAL PRODUCTS Pty Ltd

Tel : 1300 138 578
Fax: 1300 138 612

ABN 20 086 942 421

www.majacmedical.com.au

Email: majac@majacmedical.com.au