

ERKA Sphygmomanometer

Product

Description

EK223



ERKA `KOBOLD` Aneroid Sphygmomanometer Alloy Case Latex free

Adult - colours. Cat No. EK223BL = Blue, G = Blue, R = Red.

Code	Description
EK223	ERKA 'KOBOLD' Aneroid Sphygmomanometer Alloy Case, Latex Free, Colours

Code	Description
EK223B	ERKA 'SWITCH' Sphygmomanometer with Child, Infant & Newborn cuffs

Our features

- Precision manometer, anodised aluminium, Scale diameter 48mm with sensitive air-reduction valve, single tube model.
- Scale: 300mm Hg
2mm scale divisions
- New technology with air intake via the manometer capsule in the specially tempered copper-beryllium membrane.
- Black, kind to skin velcro cuff made of 100% cotton suitable for washing or dry cleaning.

Your convenience

- Elegant timeless form.
- Clear, optically large display on scale.
- Compact and robust design intended for frequent use in clinics or consulting rooms.
- Sensitively adjustable air pressure relief valve, high quality chrome finish.
- Take your pick, the models are available in black, red, blue, gold and anodised silver.
- Efficient instrument in ERKA quality for measuring blood pressure.

Kobald Sphygmomanometer with black velcro cuff single-tube model, Anodised, available in different colours.

Code no:

Silver: EK223S
Gold: EK223G
Blue: EK223BL
Red: EK223R
Black: EK223BLA

Pure Hygiene: only ERKA cuffs can be cleaned with most clinical detergents and disinfectants (e.g. Clindet & Clinicol) The sturdy materials guarantee maximum durability, and are all dermatologically tested.

EK223V010



ERKA `KOBOLD` Aneroid Sphygmomanometer

Adult, Alloy Case, Non Magnetic

Code	Description
EK223V010	ERKA 'KOBOLD' Aneroid Sphygmomanometer Adult, Alloy Case, Non Magnetic

AVAILABLE FROM:



MAJAC MEDICAL PRODUCTS Pty Ltd

Tel : 1300 138 578
Fax: 1300 138 612

ABN 20 086 942 421

www.majacmedical.com.au

Email: sales@majacmedical.com.au

ES293



ERKA - SWITCH - Sphygmomanometer - Latex free

Description

- * Comfortable handling thanks to its flat construction,
- * Outstanding ergonomics and simple adaptation for left- and right-handed use,
- * Coloured labelling and individual personalization,
- * Replaceable manometer housing,
- * Reusable handhold,
- * Calibrated cuffs, latex free
- * The service is carried out quickly
- * The technology saves both time and money
- * No depletion in the number of operational gauges available
- * No complicated logistical transport problems
- * No time lost through sending equipment
- * Minimal expenditure on documentation
- * The station remains flexible

Code	Description
ES293	ERKA 'SWITCH' Sphygmomanometer Adult, Latex Free



The brass valve screw works smoothly and is adjustable to exact settings. User-friendly valve cap for easy use, even with gloves.



The sphygmomanometer is divided into two parts. The housing can be replaced very simply. The clip fastening, which holds the two parts together, serves at the same time as the display window for the personal labelling of the gauge with colours and script.



This inflation ball has been newly designed, with the attachment neck-inserted so as to give the ball an increased volume and thereby facilitate comfortable pumping action.

By using the plastic bladder the calibrated cuff is very flexible and comfortable. The one tube technique permits production of highly comfortable cuffs, even for the smallest dimensions.

Any necessary service work on the ERKA SWITCH (e.g. replacement of the manometer housing) is carried out on site together with the MTK (Medical Technical Control). No costs due to equipment failure are incurred, the handhold is retained, so the equipment offers excellent value and economizes on resources.

The ERKA SWITCH is manufactured from shock-resistant polycarbonate. It is both elegant and robust, and has passed the shock test according to AAMI SPg. The sight glass is chemically stable and unbreakable.



A simple adjustment converts from right to left handed model



Only 'Made in Germany' is made in Germany... today and in the future!

AVAILABLE FROM:



MAJAC MEDICAL PRODUCTS Pty Ltd

Tel : 1300 138 578
Fax: 1300 138 612

ABN 20 086 942 421

www.majacmedical.com.au

Email: sales@majacmedical.com.au

EV288ML



ERKA - VARIO - Aneroid Sphyg - Light weight trolley - Latex Free

Vario basic model with cuff basket at the front and spiral tube, for universal use + Trolley.

- * high quality precision manometer,
 - * diameter of scale 150 mm, with mirrored exterior ring,
 - * scale inclined by 12°,
 - * housing made of high-quality ABS in signal white/light-grey,
 - * non-reflective manometer-glass,
 - * cuff-basket in powder-coated steel,
 - * inflation bulb with sensitive adjustable screw reduction valve,
- The Vario basic model combines with simple accessories for universal use:

- * Desk model, required accessory: pedestal,
- * Wall model with swivel motion required accessory: wall bracket,
- * Rail-mounted model with swivel motion, required accessory: universal clamp height 15-50 mm, * Stand model on 5 castors required accessory: stand.
- * choice of 8 different cuff-sizes, velcro cuff:
 - cotton poplin black or green,
 - PUR coated fabric, washable,
- colours: green, red or blue,
- * all parts replaceable,
- * 36 months guarantee against faulty materials and workmanship

Code	Description
EV288ML	ERKA 'VARIO' mobile Aneroid Sphygmomanometer, light weight mobile stand

EV288WA



ERKA `VARIO` Wall model - Aneroid Adult - Latex free

ERKA - VARIO - Aneroid Sphyg - Light weight trolley - Latex Free Vario basic model with cuff basket at the front and spiral tube, for universal use + Trolley.

- * high quality precision manometer,
 - * diameter of scale 150 mm, with mirrored exterior ring,
 - * scale inclined by 12°,
 - * housing made of high-quality ABS in signal white/light-grey,
 - * non-reflective manometer-glass,
 - * cuff-basket in powder-coated steel,
 - * inflation bulb with sensitive adjustable screw reduction valve,
- The Vario basic model combines with simple accessories for universal use:

- * Desk model, required accessory: pedestal,
- * Wall model with swivel motion required accessory: wall bracket,
- * Rail-mounted model with swivel motion, required accessory: universal clamp height 15-50 mm, * Stand model on 5 castors required accessory: stand.
- * choice of 8 different cuff-sizes, velcro cuff:
 - cotton poplin black or green,
 - PUR coated fabric, washable,
- colours: green, red or blue,
- * all parts replaceable,
- * 36 months guarantee against faulty materials and workmanship

Code	Description
EV288WA	ERKA 'VARIO' wall model, Aneroid Adult, Latex Free

AVAILABLE FROM:



MAJAC MEDICAL PRODUCTS Pty Ltd

Tel : 1300 138 578
Fax: 1300 138 612

ABN 20 086 942 421

www.majacmedical.com.au

Email: sales@majacmedical.com.au

EE176DMM



ERKA Erkameter 3000 - Desk Model

Mercurial Sphyg Adult - Unbreakable ABS Case

Description:

- * portable desk model in a high-quality plastic case,
- * unbreakable capillary system,
- * square section capillary to read out any measured value without parallax error,
- * more than 40 percent mercury saved for ecological reasons,
- * automatic mercury safety catch,
- * modern ergonomic design,
- * reduction valve for sensitive adjustment, high-quality chromium-plated. * all parts replaceable,
- * 36 months guarantee against faulty materials and workmanship.

Code	Description
EE176DMM	ERKA desk model 3000 – Mercurial Sphygmomanometer Adult, Unbreakable ABS Case

EV285



ERKA Vario Desk model with basket at the back

Description

- * high quality precision manometer,
- * diameter of scale 150 mm, with mirrored exterior ring,
- * scale inclined by 12°,
- * housing made of high-quality ABS in signal white/light-grey,
- * non-reflective manometer-glass,
- * cuff-basket in powder-coated steel,
- * inflation bulb with sensitive adjustable screw reduction valve, velcro cuff:
 - cotton poplin black or green,
 - PUR coated fabric, washable,
- colours: green, red or blue,
- * all parts replaceable, * 36 months guarantee against faulty materials and workmanship

Code	Description
EV285	ERKA 'VARIO' Desk Aneroid, Adult, Latex Free

EVRM295



Rail Mount for VARIO Sphygmomanometer

Code	Description
EVRM295	Rail Mount for VARIO Sphygmomanometer

AVAILABLE FROM:



MAJAC MEDICAL PRODUCTS Pty Ltd

Tel : 1300 138 578

Fax: 1300 138 612

ABN 20 086 942 421

www.majacmedical.com.au

Email: sales@majacmedical.com.au