

# ASKIR SURGICAL ASPIRATOR

**ASKIR 230-12 BR Code: SPASKIRBR**  
**Suction Pump 240/12 Volt with battery**

The new **ASKIR 230-12 BR** is a device working with 230V ~ AC 12V DC mixed power supply. Easily transportable from one location to another - particularly suitable for tracheotomized patients, minor surgical applications and post-operative therapy at home. To be used for the nasal, oral & tracheal aspiration of body liquids (mucus or catarrh) in adults and children. Designed for almost continual use, thanks to an electronic system for managing the electrical supply and monitoring the functioning of the internal battery - enabling it to be used at a voltage of 230 V even when the battery is being recharged. A series of LED on the control panel signal whether the machine is operating, as well as the status of the batteries.



Made of highly heat resistant, electrically insulated plastic in conformity with the latest European safety standard. The unit is supplied with a complete polycarbonate autoclavable jar with overflow valve and it is equipped with aspiration regulator and vacuum indicator.

## Technical Characteristics:

<b>Voltage:</b>	230V/50Hz 12 Volt DC Car Adaptor Internal Battery	<b>Weight:</b>	3.6kg
<b>Power Consumption:</b>	50VA	<b>Dimensions:</b>	350(w) x 210(h) x 180mm(d)
<b>Fuse:</b>	1 x 1.6 A 250V	<b>Standard Accessories:</b>	Tubes Set 6 x 10mm Probe Connector Aspiration Probe CH20 Antibacterial Filter Aspiration Jar 1000cc
<b>ISO 10079-1:</b>	High Vacuum/Low Flow		
<b>Maximum suction aspiration:</b>	-0.08 Bar		
<b>Maximum suction flow:</b>	16 l/min		
<b>Battery Life:</b>	45-50 minutes		
<b>Battery Charge Time:</b>	120-150 minutes		

AVAILABLE FROM:



**MAJAC MEDICAL PRODUCTS Pty Ltd**

ABN 20 086 942 421

Tel : 1300 138 578

[www.majacmedical.com.au](http://www.majacmedical.com.au)

Fax: 1300 138 612

Email: [majac@optusnet.com.au](mailto:majac@optusnet.com.au)

# ASKIR SURGICAL ASPIRATOR



**ASKIR 30** Code: **SPASKIR30**  
**Suction Pump 30 Litre 240 Volt**

The new **ASKIR 30** is a surgical aspirator working 230V/50Hz network electricity, to be used for nasal, oral & tracheal aspiration of body liquids (mucus, catarrh or blood) in adults or children. An easily portable machine designed for non continuous use of 120 min ON / 60 min OFF. Ideal hospital use for tracheotomized patients, minor surgical applications and post-operative therapy at home.

**Weight:** 3.6kg  
**Dimensions:** 350(w) x 210(h) x 180mm(d)  
**Standard Accessories:** Tubes Set 6x10mm  
Probe Connector  
Aspiration Probe CH20  
Antibacterial Filter  
Aspiration Jar 1000cc



**ASKIR 20** Code: **SPASKIR20**  
**Suction Pump 16 Litre 240 Volt**

The new **ASKIR 20** is a surgical aspirator working 230V/50Hz network electricity, to be used for nasal, oral & tracheal aspiration of body liquids (mucus, catarrh or blood) in adults or children. An easily portable machine designed for non continuous use of 30 min ON / 30 min OFF. Ideal hospital use for tracheotomized patients, minor surgical applications and post-operative therapy at home.

**Weight:** 2.5kg  
**Dimensions:** 350(w) x 210(h) x 180mm(d)  
**Standard Accessories:** Tubes Set 6x10mm  
Probe Connector  
Aspiration Probe CH20  
Antibacterial Filter  
Aspiration Jar 1000cc

AVAILABLE FROM:



**MAJAC MEDICAL PRODUCTS Pty Ltd**

ABN 20 086 942 421

[www.majacmedical.com.au](http://www.majacmedical.com.au)

Email: [majac@optusnet.com.au](mailto:majac@optusnet.com.au)

Tel : 1300 138 578

Fax: 1300 138 612