

# No Touch






FEVER MEASUREMENT WITHOUT CONTACT



**IDEAL FOR UNCOOPERATIVE PATIENTS & CHILDREN**  
(NO NEED TO WAKE OR TOUCH PATIENT)



## Features and Specifications

-  Automatic and comfortable measurement by scanning the temple area from a distance of up to 2.5cm.
-  Measuring within seconds.
-  8 memory recall.
-  Accuracy according to ASTM standard.
-  No Touch reduces transmission of infections.



**MAJAC MEDICAL PRODUCTS** Pty. Ltd.

Sales inquiries TEL: 1300 138 578 FAX: 1300 138 612  
Business Office TEL: +61 7 3265 6355 FAX: +61 7 3265 2729  
E-Mail: [sales@majacmedical.com.au](mailto:sales@majacmedical.com.au)  
[www.majacmedical.com.au](http://www.majacmedical.com.au)

*Your Complete Range of Cleaning & Infection Control Products*

# No Touch Thermometer

Fever measurement  
without contact

## Temporal artery thermometry



**Hygienic, safe and economic.**



**User and patient friendly**



**Suitable from paediatrics to geriatrics**



**Ideal for uncooperative patients and children,  
no need to stand still, wake-up or touch.**



**Currently servicing hospitals, aged care facilities  
and many GP's**



## How does it work ?

Using advanced technology the “**No Touch**” measures the temperature of the infra red heat given off by the temple skin and converts this temperature into an oral equivalent temperature.



**MAJAC MEDICAL PRODUCTS** Pty. Ltd.

Sales Inquiries TEL: 1300 138 578 FAX: 1300 138 612

Business Office TEL: +61 7 3265 6355 FAX: +61 7 3265 2729

E-Mail: [sales@majacmedical.com.au](mailto:sales@majacmedical.com.au)

[www.majacmedical.com.au](http://www.majacmedical.com.au)

*Your Complete Range of Cleaning & Infection Control Products*