

DEFIBRILLATORS

CU-SP1 (SP1 - IPAD AED)

Easy-to-use Semi-Automated External Defibrillator (AED) that is small, light, and portable

KEY FEATURES:

- Designed to be used by minimally trained individuals
- Ambient noise detector measures level of background noise and adjusts the volume of the voice prompts accordingly
- Integrated pad storage: Pre-connected pads are stored in a clear compartment on the underside of the unit
- Visual indicator shows at a glance battery life, unit status and pad status

OPTIONS AND ACCESSORIES:

CAT. No.	DESCRIPTION
CU-SP1-A-BAG-3010	SPI Carrying case
CU-A1007S	Disposable Defib pads (Adult)
CU-A1102S	Disposable Defib pads (Child)
CU-SPI-A08	Wall cabinet
CU-SPI-A12	Wall bracket
CU-SA1103BB	Disposable battery pack
CU-SPI-A10-SD	SD card
CU-SPT	Trainer model

- Easily switch from adult to child mode without changing pads
- Auto volume adjustment
- CPR detection
- Daily, weekly, monthly self checks



Excludes Battery



CU-SP1

SP1-IPAD SP1

CAT. No.	DESCRIPTION
CU-SP1	iPAD SP1 - Semi Automated External Defibrillator AED)



Wall Cabinet



Carry Case
CU-SP1-A-BAG-3010



Wall Bracket



Smart Pad Storage



Smart CPR Detection



Majac Medical Products Pty Ltd
 Ph: +617 3265 6355
 Fax: +617 3865 2729
 sales@majacmedical.com.au
 www.majacmedical.com.au



FS 540017

Also Available From: